



URGENT MEDICAL DEVICE CORRECTION

Heart Failure Division
Thoratec Switzerland GmbH
Technoparkstrasse 1
CH-8005 Zürich Switzerland

2nd Generation CentriMag Primary Console
(Model Number: 201-90411, L201-90411)
(UDI: 07640135140702, 07640135140702)

May 2024

Dear Valued Customer,

Abbott is notifying users of a potential issue identified through our internal process with the **CentriMag Console (2nd Generation CentriMag Primary Console)** which is a part of the CentriMag™ Circulatory Support System. Internal testing showed that due to a component change, the CentriMag Console may not conform to the IEC 61000-4-5 standard. If subjected to a power surge above 1.8kV and up to 2.0 kV, the CentriMag Console may shut down completely without alarm. A power surge may be caused by, including, but not limited to the following: lightning strike, unregulated power switching, or faulty equipment. Pumping would cease and the patient would have to be switched to a backup.

There have been no reported complaints for distributed products and no reports of patient harm or adverse events due to this issue. Only surge voltage above 1.8 kV may shut down the CentriMag console without alarm. The ability of the user to transfer to a backup console to restore system function is not impacted. The affected product is not being removed from the field and does not need to be returned.

This letter contains important information and recommendations on how to use CentriMag Console if the system shuts down without an alarm.

Abbott discovered that certain CentriMag Console manufactured since March 21, 2016 are impacted by this issue. Please review Appendix A for the specific model and serial numbers within the scope of this notification.

Impact and Associated Risks

When affected CentriMag Console is exposed to surge voltage above 1.8 kV and up to 2.0 kV, if the Console shuts down without alarming, resulting impact can be delay in hemodynamic support while support is transferred to the backup system, or potentially thromboembolism if patient support is not transferred in a timely manner.

Recommendation

If the CentriMag Console shuts down without alarms or alerts present, follow Abbott Instruction for Use (IFU) for *Console and Motor malfunction* as described in the Section “Emergency / Troubleshooting”, and switch to a backup CentriMag Console which must be prepared as required backup components in the immediate vicinity.

Note: When the backup CentriMag console is put into use, the previously used console can be restarted and set up as a backup console.

Per the IFU and current practice, a full CentriMag Console backup system is always required to be available in the immediate vicinity of each patient whenever the CentriMag system is used.



Abbott

URGENT MEDICAL DEVICE CORRECTION

Heart Failure Division
Thoratec Switzerland GmbH
Technoparkstrasse 1
CH-8005 Zürich Switzerland

2nd Generation CentriMag Primary Console
(Model Number: 201-90411, L201-90411)
(UDI: 07640135140702, 07640135140702)

Abbott performed additional testing and confirmed that the ability to pump and alarm is not affected in backup systems. As stated in the IFU these backup systems should have a battery charge sufficient for at least one hour of operation, be connected to wall power, and be powered off.

Abbott is in the process of implementing an additional in-process test by Q3 of 2024, which will ensure that all systems which will be distributed are confirmed to withstand a surge up to 2.0 kV.

Please complete the acknowledgement form included with this letter and return to Abbott. Please distribute this notice to those who need to be aware within your institution and to any organization where potentially affected devices may have been transferred in the event that devices were transferred elsewhere.

Abbott has notified applicable regulatory agencies about this issue. Adverse reactions or quality problems experienced with the use of this product may be reported to the FDA's MedWatch Adverse Event Reporting program either online, by regular mail or by fax. To submit your report:

- Complete the voluntary Form FDA 3500 online
- Download form from [FDA.gov](https://www.fda.gov) or call 1-800-332-1088 to request a reporting form, then complete and return to the address on the pre-addressed form or submit by fax to 1-800-FDA-0178.

Abbott is committed to providing the highest quality products and support, and we thank you for your partnership in assisting us with this process. Should you have any questions about this communication, please contact Abbott Customer Service at 1 (800) 456-1477 or your local Abbott representative.

Sincerely,

Dan Lack
Senior Manager Quality
Quality Management Representative
Abbott Heart Failure

Appendix A: Specific Model and Serial Numbers in Scope of Notification

CentriMag Console Serial Number Range: L04549-0006, 19234128 to 19359336

Model Number	Serial Number
201-30300	L04549-0006
201-90411	19234128
201-90411	19250285
201-90411	19260491
201-90411	19260493
201-90411	19260849
201-90411	19260850
201-90411	19260869
201-90411	19260870
201-90411	19260871
201-90411	19260872
201-90411	19260873
201-90411	19260874
201-90411	19260875
201-90411	19260876
201-90411	19260877
201-90411	19260928
201-90411	19261079
201-90411	19261081
201-90411	19261082
201-90411	19261083
201-90411	19261084
201-90411	19261085
201-90411	19261087
201-90411	19261088
201-90411	19261089
201-90411	19261090
201-90411	19261091
201-90411	19261092
201-90411	19271685
201-90411	19271686
201-90411	19271687
201-90411	19271688
201-90411	19271689
201-90411	19272717
201-90411	19272718
201-90411	19272719
201-90411	19276573
201-90411	19276574
201-90411	19277457
201-90411	19277458
201-90411	19277459
201-90411	19277460
201-90411	19277461
201-90411	19279320
201-90411	19279321

Model Number	Serial Number
201-90411	19279322
201-90411	19280692
201-90411	19280693
201-90411	19280694
201-90411	19280696
201-90411	19280697
201-90411	19280698
201-90411	19280699
201-90411	19280700
201-90411	19280701
201-90411	19280703
201-90411	19280704
201-90411	19280705
201-90411	19280706
201-90411	19284718
201-90411	19284719
201-90411	19284720
201-90411	19284721
201-90411	19284722
201-90411	19284723
201-90411	19284724
201-90411	19284725
201-90411	19284726
201-90411	19284727
201-90411	19284728
201-90411	19284729
201-90411	19284730
201-90411	19284731
201-90411	19284732
201-90411	19284733
201-90411	19284734
201-90411	19284735
201-90411	19284736
201-90411	19284737
201-90411	19284738
201-90411	19284739
201-90411	19284741
201-90411	19284742
201-90411	19284743
201-90411	19284744
201-90411	19284745
201-90411	19284746
201-90411	19303709
201-90411	19303710
201-90411	19303711
201-90411	19303712

Model Number	Serial Number
201-90411	19303713
201-90411	19303714
201-90411	19303715
201-90411	19303716
201-90411	19303717
201-90411	19303718
201-90411	19303719
201-90411	19303720
201-90411	19303721
201-90411	19303722
201-90411	19303723
201-90411	19303725
201-90411	19303726
201-90411	19303727
201-90411	19303728
201-90411	19303729
201-90411	19303730
201-90411	19303731
201-90411	19303732
201-90411	19303733
201-90411	19303734
201-90411	19303735
201-90411	19303736
201-90411	19303737
201-90411	19303738
201-90411	19303739
201-90411	19303740
201-90411	19303741
201-90411	19303742
201-90411	19303743
201-90411	19303746
201-90411	19303747
201-90411	19336185
201-90411	19336188
201-90411	19343549
201-90411	19343550
201-90411	19343551
201-90411	19343552
201-90411	19355272
201-90411	19355273
201-90411	19355274
201-90411	19355276
201-90411	19359270
201-90411	19359281
201-90411	19359302
201-90411	19359336

Appendix A: Specific Model and Serial Numbers in Scope of Notification

CentriMag Console Serial Number Range: 19359337 to 20076171

Model Number	Serial Number
201-90411	19359337
201-90411	19359338
201-90411	19361777
201-90411	19361778
201-90411	19361779
201-90411	19361781
201-90411	19361783
201-90411	19361784
201-90411	19361785
201-90411	19366820
201-90411	19366821
201-90411	19366829
201-90411	19366830
201-90411	19366831
201-90411	19366932
201-90411	19366936
201-90411	19366939
201-90411	19366940
201-90411	19366941
201-90411	19366942
201-90411	19366943
201-90411	19369083
201-90411	19369084
201-90411	19369086
201-90411	19372718
201-90411	19379714
201-90411	19379715
201-90411	19379730
201-90411	19383630
201-90411	19383635
201-90411	19383636
201-90411	19383637
201-90411	19400415
201-90411	19400416
201-90411	19400417
201-90411	19400418
201-90411	19400423
201-90411	19404180
201-90411	19404181
201-90411	19407985
201-90411	19411156
201-90411	19421179
201-90411	19429363
201-90411	19430948
201-90411	19430950
201-90411	19430952

Model Number	Serial Number
201-90411	19430953
201-90411	19434048
201-90411	19434049
201-90411	19434053
201-90411	19434054
201-90411	19434056
201-90411	19450851
201-90411	19450852
201-90411	19450853
201-90411	19450854
201-90411	19450855
201-90411	19450856
201-90411	19450857
201-90411	19450858
201-90411	19450861
201-90411	19450862
201-90411	19474334
201-90411	19507820
201-90411	19535368
201-90411	19557751
201-90411	19604015
201-90411	19604018
201-90411	19605456
201-90411	19605457
201-90411	19923332
201-90411	19923534
201-90411	19924934
201-90411	19924935
201-90411	19926286
201-90411	19926287
201-90411	19926289
201-90411	19933471
201-90411	19942145
201-90411	19942150
201-90411	19944575
201-90411	19944577
201-90411	19944578
201-90411	19944580
201-90411	19950806
201-90411	19950809
201-90411	19955809
201-90411	19955810
201-90411	19955811
201-90411	19955814
201-90411	19960230
201-90411	19960231

Model Number	Serial Number
201-90411	19960232
201-90411	19960233
201-90411	19960234
201-90411	19969063
201-90411	19969065
201-90411	19981478
201-90411	19981955
201-90411	19981956
201-90411	19981958
201-90411	19981959
201-90411	19989560
201-90411	19989561
201-90411	19990920
201-90411	19990921
201-90411	19990932
201-90411	19990933
201-90411	19999967
201-90411	20001637
201-90411	20001691
201-90411	20001692
201-90411	20001693
201-90411	20035494
201-90411	20035495
201-90411	20035496
201-90411	20035497
201-90411	20035498
201-90411	20052389
201-90411	20052390
201-90411	20059198
201-90411	20059643
201-90411	20060491
201-90411	20060494
201-90411	20064025
201-90411	20064689
201-90411	20064690
201-90411	20064691
201-90411	20070964
201-90411	20070965
201-90411	20070966
201-90411	20070968
201-90411	20073850
201-90411	20073851
201-90411	20073852
201-90411	20073853
201-90411	20076100
201-90411	20076171

Appendix A: Specific Model and Serial Numbers in Scope of Notification

CentriMag Console Serial Number Range: 20077699 to 20190463, L03849-0004 to L04696-0007

Model Number	Serial Number
201-90411	20077699
201-90411	20077700
201-90411	20077701
201-90411	20078986
201-90411	20079883
201-90411	20084412
201-90411	20084414
201-90411	20087457
201-90411	20087458
201-90411	20087459
201-90411	20092097
201-90411	20092098
201-90411	20093315
201-90411	20093317
201-90411	20093763
201-90411	20093764
201-90411	20108705
201-90411	20111289
201-90411	20114030
201-90411	20114031
201-90411	20114032
201-90411	20114033
201-90411	20114034
201-90411	20116579
201-90411	20116581
201-90411	20118497
201-90411	20119744
201-90411	20119745
201-90411	20119746
201-90411	20119747
201-90411	20124032
201-90411	20126798
201-90411	20126799
201-90411	20126800
201-90411	20134795
201-90411	20134796
201-90411	20134797
201-90411	20137425
201-90411	20137428
201-90411	20137429
201-90411	20139856
201-90411	20139857
201-90411	20139858
201-90411	20139859
201-90411	20141944
201-90411	20141945

Model Number	Serial Number
201-90411	20144706
201-90411	20144708
201-90411	20144709
201-90411	20154839
201-90411	20157604
201-90411	20160894
201-90411	20160895
201-90411	20164796
201-90411	20164797
201-90411	20164799
201-90411	20164801
201-90411	20164804
201-90411	20164805
201-90411	20164806
201-90411	20164807
201-90411	20181207
201-90411	20181209
201-90411	20184744
201-90411	20188003
201-90411	20188360
201-90411	20190463
201-90411	L03849-0004
201-90411	L03849-0009
201-90411	L03897-0001
201-90411	L03897-0002
201-90411	L04008-0001
201-90411	L04008-0002
201-90411	L04209-0007
201-90411	L04209-0008
201-90411	L04342-0001
201-90411	L04342-0004
201-90411	L04342-0006
201-90411	L04342-0008
201-90411	L04350-0003
201-90411	L04350-0005
201-90411	L04350-0009
201-90411	L04350-0010
201-90411	L04380-0001
201-90411	L04380-0002
201-90411	L04380-0003
201-90411	L04380-0004
201-90411	L04380-0005
201-90411	L04380-0006
201-90411	L04380-0007
201-90411	L04380-0008
201-90411	L04380-0009

Model Number	Serial Number
201-90411	L04380-0010
201-90411	L04406-0001
201-90411	L04406-0004
201-90411	L04437-0001
201-90411	L04437-0002
201-90411	L04437-0003
201-90411	L04437-0004
201-90411	L04437-0005
201-90411	L04437-0006
201-90411	L04437-0007
201-90411	L04437-0008
201-90411	L04437-0009
201-90411	L04503-0001
201-90411	L04503-0002
201-90411	L04512-0006
201-90411	L04549-0005
201-90411	L04549-0007
201-90411	L04549-0008
201-90411	L04549-0009
201-90411	L04549-0010
201-90411	L04556-0001
201-90411	L04556-0002
201-90411	L04556-0003
201-90411	L04594-0002
201-90411	L04594-0003
201-90411	L04594-0004
201-90411	L04594-0005
201-90411	L04594-0006
201-90411	L04594-0007
201-90411	L04648-0001
201-90411	L04648-0002
201-90411	L04648-0003
201-90411	L04685-0002
201-90411	L04685-0003
201-90411	L04685-0005
201-90411	L04685-0006
201-90411	L04685-0007
201-90411	L04685-0008
201-90411	L04685-0009
201-90411	L04685-0010
201-90411	L04696-0001
201-90411	L04696-0002
201-90411	L04696-0003
201-90411	L04696-0004
201-90411	L04696-0006
201-90411	L04696-0007

Appendix A: Specific Model and Serial Numbers in Scope of Notification

CentriMag Console Serial Number Range: L04696-0008 to L05254-0001

Model Number	Serial Number
201-90411	L04696-0008
201-90411	L04696-0009
201-90411	L04696-0010
201-90411	L04729-0001
201-90411	L04729-0002
201-90411	L04729-0003
201-90411	L04729-0004
201-90411	L04729-0005
201-90411	L04729-0007
201-90411	L04729-0008
201-90411	L04729-0009
201-90411	L04729-0010
201-90411	L04734-0003
201-90411	L04734-0004
201-90411	L04734-0005
201-90411	L04734-0008
201-90411	L04734-0009
201-90411	L04734-0010
201-90411	L04742-0001
201-90411	L04742-0003
201-90411	L04742-0004
201-90411	L04742-0005
201-90411	L04742-0007
201-90411	L04753-0006
201-90411	L04753-0007
201-90411	L04753-0008
201-90411	L04753-0009
201-90411	L04753-0010
201-90411	L04766-0001
201-90411	L04766-0002
201-90411	L04766-0003
201-90411	L04766-0004
201-90411	L04766-0005
201-90411	L04766-0006
201-90411	L04766-0007
201-90411	L04766-0008
201-90411	L04766-0009
201-90411	L04766-0010
201-90411	L04776-0001
201-90411	L04776-0002
201-90411	L04776-0003
201-90411	L04776-0005
201-90411	L04800-0001
201-90411	L04800-0002
201-90411	L04800-0003
201-90411	L04800-0004

Model Number	Serial Number
201-90411	L04800-0005
201-90411	L04800-0006
201-90411	L04800-0007
201-90411	L04800-0008
201-90411	L04800-0009
201-90411	L04800-0010
201-90411	L04852-0001
201-90411	L04852-0002
201-90411	L04852-0003
201-90411	L04852-0004
201-90411	L04852-0005
201-90411	L04852-0006
201-90411	L04852-0007
201-90411	L04852-0008
201-90411	L04852-0009
201-90411	L04852-0010
201-90411	L04862-0001
201-90411	L04862-0002
201-90411	L04862-0003
201-90411	L04862-0004
201-90411	L04862-0005
201-90411	L04862-0006
201-90411	L04862-0007
201-90411	L04862-0008
201-90411	L04862-0009
201-90411	L04862-0010
201-90411	L04884-0002
201-90411	L04884-0005
201-90411	L04884-0008
201-90411	L04884-0009
201-90411	L04884-0010
201-90411	L04900-0001
201-90411	L04900-0002
201-90411	L04900-0003
201-90411	L04900-0004
201-90411	L04900-0005
201-90411	L04923-0001
201-90411	L04923-0002
201-90411	L04923-0007
201-90411	L04923-0008
201-90411	L04923-0010
201-90411	L04942-0008
201-90411	L04942-0009
201-90411	L04942-0010
201-90411	L04966-0005
201-90411	L04966-0008

Model Number	Serial Number
201-90411	L04966-0010
201-90411	L04982-0001
201-90411	L04982-0002
201-90411	L05083-0008
201-90411	L05083-0009
201-90411	L05083-0010
201-90411	L05096-0001
201-90411	L05116-0001
201-90411	L05116-0002
201-90411	L05116-0003
201-90411	L05116-0004
201-90411	L05116-0005
201-90411	L05116-0010
201-90411	L05129-0003
201-90411	L05129-0004
201-90411	L05129-0005
201-90411	L05129-0006
201-90411	L05150-0001
201-90411	L05150-0002
201-90411	L05150-0005
201-90411	L05168-0001
201-90411	L05168-0002
201-90411	L05168-0003
201-90411	L05168-0004
201-90411	L05168-0005
201-90411	L05168-0009
201-90411	L05202-0001
201-90411	L05202-0002
201-90411	L05202-0003
201-90411	L05202-0004
201-90411	L05202-0007
201-90411	L05212-0003
201-90411	L05212-0004
201-90411	L05212-0005
201-90411	L05212-0006
201-90411	L05212-0007
201-90411	L05229-0002
201-90411	L05229-0003
201-90411	L05229-0004
201-90411	L05229-0005
201-90411	L05229-0006
201-90411	L05241-0004
201-90411	L05241-0006
201-90411	L05241-0007
201-90411	L05241-0010
201-90411	L05254-0001

Appendix A: Specific Model and Serial Numbers in Scope of Notification

CentriMag Console Serial Number Range: L05254-0002 to L06066-0003

Model Number	Serial Number
201-90411	L05254-0002
201-90411	L05254-0003
201-90411	L05254-0004
201-90411	L05308-0001
201-90411	L05326-0004
201-90411	L05326-0005
201-90411	L05326-0006
201-90411	L05326-0007
201-90411	L05326-0009
201-90411	L05326-0010
201-90411	L05416-0006
201-90411	L05416-0007
201-90411	L05416-0008
201-90411	L05416-0009
201-90411	L05416-0010
201-90411	L05427-0001
201-90411	L05427-0002
201-90411	L05427-0003
201-90411	L05427-0004
201-90411	L05427-0005
201-90411	L05427-0006
201-90411	L05427-0007
201-90411	L05437-0008
201-90411	L05437-0009
201-90411	L05437-0010
201-90411	L05439-0001
201-90411	L05439-0002
201-90411	L05439-0003
201-90411	L05451-0001
201-90411	L05451-0002
201-90411	L05451-0003
201-90411	L05451-0004
201-90411	L05451-0005
201-90411	L05451-0006
201-90411	L05451-0007
201-90411	L05451-0008
201-90411	L05451-0010
201-90411	L05569-0003
201-90411	L05569-0004
201-90411	L05569-0005
201-90411	L05569-0006
201-90411	L05569-0008
201-90411	L05569-0009
201-90411	L05575-0001
201-90411	L05575-0002
201-90411	L05575-0003

Model Number	Serial Number
201-90411	L05575-0004
201-90411	L05575-0005
201-90411	L05575-0006
201-90411	L05575-0007
201-90411	L05575-0010
201-90411	L05582-0001
201-90411	L05582-0002
201-90411	L05582-0003
201-90411	L05582-0004
201-90411	L05582-0005
201-90411	L05582-0007
201-90411	L05582-0008
201-90411	L05582-0010
201-90411	L05618-0001
201-90411	L05618-0003
201-90411	L05639-0004
201-90411	L05639-0005
201-90411	L05639-0006
201-90411	L05639-0007
201-90411	L05639-0008
201-90411	L05639-0009
201-90411	L05639-0010
201-90411	L05654-0001
201-90411	L05654-0004
201-90411	L05708-0001
201-90411	L05708-0002
201-90411	L05708-0003
201-90411	L05708-0004
201-90411	L05708-0005
201-90411	L05708-0006
201-90411	L05708-0007
201-90411	L05708-0008
201-90411	L05708-0009
201-90411	L05708-0010
201-90411	L05724-0001
201-90411	L05724-0002
201-90411	L05724-0006
201-90411	L05724-0007
201-90411	L05724-0009
201-90411	L05724-0010
201-90411	L05729-0002
201-90411	L05729-0003
201-90411	L05729-0004
201-90411	L05729-0007
201-90411	L05729-0008
201-90411	L05729-0009

Model Number	Serial Number
201-90411	L05729-0010
201-90411	L05779-0002
201-90411	L05779-0007
201-90411	L05779-0008
201-90411	L05837-0001
201-90411	L05837-0002
201-90411	L05837-0003
201-90411	L05837-0004
201-90411	L05837-0007
201-90411	L05837-0008
201-90411	L05837-0010
201-90411	L05851-0002
201-90411	L05851-0003
201-90411	L05851-0004
201-90411	L05851-0005
201-90411	L05851-0006
201-90411	L05851-0007
201-90411	L05851-0008
201-90411	L05851-0009
201-90411	L05851-0010
201-90411	L05877-0004
201-90411	L05877-0005
201-90411	L05877-0006
201-90411	L05877-0007
201-90411	L05917-0002
201-90411	L05917-0004
201-90411	L05917-0005
201-90411	L05917-0007
201-90411	L05963-0010
201-90411	L05966-0001
201-90411	L05966-0002
201-90411	L05966-0003
201-90411	L06060-0008
201-90411	L06060-0009
201-90411	L06060-0010
201-90411	L06064-0001
201-90411	L06064-0002
201-90411	L06064-0003
201-90411	L06064-0004
201-90411	L06064-0005
201-90411	L06064-0006
201-90411	L06064-0009
201-90411	L06064-0010
201-90411	L06066-0001
201-90411	L06066-0002
201-90411	L06066-0003

Appendix A: Specific Model and Serial Numbers in Scope of Notification

CentriMag Console Serial Number Range: L06066-0004 to L06401-0001

Model Number	Serial Number
201-90411	L06066-0004
201-90411	L06066-0005
201-90411	L06066-0006
201-90411	L06066-0007
201-90411	L06066-0008
201-90411	L06299-0008
201-90411	L06299-0009
201-90411	L06302-0001
201-90411	L06302-0002
201-90411	L06302-0003
201-90411	L06302-0004
201-90411	L06302-0005
201-90411	L06302-0006
201-90411	L06302-0007
201-90411	L06302-0009
201-90411	L06302-0010
201-90411	L06308-0001
201-90411	L06308-0002
201-90411	L06308-0003
201-90411	L06308-0004
201-90411	L06308-0006
201-90411	L06308-0007
201-90411	L06308-0008
201-90411	L06308-0009
201-90411	L06309-0001
201-90411	L06309-0002
201-90411	L06336-0006
201-90411	L06336-0007
201-90411	L06336-0008
201-90411	L06336-0009
201-90411	L06336-0010
201-90411	L06338-0005
201-90411	L06338-0006
201-90411	L06338-0010
201-90411	L06339-0001
201-90411	L06339-0002
201-90411	L06339-0003
201-90411	L06339-0004
201-90411	L06339-0005
201-90411	L06339-0007
201-90411	L06339-0008
201-90411	L06339-0009
201-90411	L06339-0010
201-90411	L06340-0001
201-90411	L06340-0003
201-90411	L06340-0004

Model Number	Serial Number
201-90411	L06340-0005
201-90411	L06340-0006
201-90411	L06340-0008
201-90411	L06340-0009
201-90411	L06340-0010
201-90411	L06355-0002
201-90411	L06355-0005
201-90411	L06355-0007
201-90411	L06355-0008
201-90411	L06355-0009
201-90411	L06355-0010
201-90411	L06356-0001
201-90411	L06356-0002
201-90411	L06356-0003
201-90411	L06356-0004
201-90411	L06356-0005
201-90411	L06356-0006
201-90411	L06356-0007
201-90411	L06356-0008
201-90411	L06356-0009
201-90411	L06356-0010
201-90411	L06359-0001
201-90411	L06359-0002
201-90411	L06359-0003
201-90411	L06359-0005
201-90411	L06359-0006
201-90411	L06359-0007
201-90411	L06359-0008
201-90411	L06359-0009
201-90411	L06359-0010
201-90411	L06364-0001
201-90411	L06364-0002
201-90411	L06364-0007
201-90411	L06364-0009
201-90411	L06364-0010
201-90411	L06379-0007
201-90411	L06379-0009
201-90411	L06380-0001
201-90411	L06380-0007
201-90411	L06380-0008
201-90411	L06380-0009
201-90411	L06380-0010
201-90411	L06381-0001
201-90411	L06381-0002
201-90411	L06381-0003
201-90411	L06381-0004

Model Number	Serial Number
201-90411	L06382-0001
201-90411	L06382-0002
201-90411	L06382-0003
201-90411	L06382-0004
201-90411	L06382-0005
201-90411	L06382-0006
201-90411	L06382-0007
201-90411	L06382-0008
201-90411	L06382-0009
201-90411	L06382-0010
201-90411	L06383-0001
201-90411	L06383-0002
201-90411	L06383-0008
201-90411	L06383-0009
201-90411	L06383-0010
201-90411	L06384-0001
201-90411	L06384-0002
201-90411	L06384-0003
201-90411	L06384-0004
201-90411	L06384-0005
201-90411	L06384-0006
201-90411	L06384-0007
201-90411	L06384-0008
201-90411	L06384-0009
201-90411	L06384-0010
201-90411	L06385-0001
201-90411	L06385-0002
201-90411	L06385-0003
201-90411	L06385-0010
201-90411	L06386-0001
201-90411	L06386-0002
201-90411	L06386-0004
201-90411	L06386-0005
201-90411	L06386-0006
201-90411	L06387-0001
201-90411	L06387-0002
201-90411	L06387-0003
201-90411	L06387-0004
201-90411	L06387-0007
201-90411	L06387-0008
201-90411	L06387-0009
201-90411	L06387-0010
201-90411	L06388-0005
201-90411	L06388-0006
201-90411	L06388-0009
201-90411	L06401-0001

Appendix A: Specific Model and Serial Numbers in Scope of Notification

CentriMag Console Serial Number Range: L06401-0006 to L06524-0010

Model Number	Serial Number
201-90411	L06401-0006
201-90411	L06401-0007
201-90411	L06401-0008
201-90411	L06401-0009
201-90411	L06401-0010
201-90411	L06402-0002
201-90411	L06403-0002
201-90411	L06403-0003
201-90411	L06403-0004
201-90411	L06403-0005
201-90411	L06403-0006
201-90411	L06403-0007
201-90411	L06403-0009
201-90411	L06403-0010
201-90411	L06408-0002
201-90411	L06408-0003
201-90411	L06408-0004
201-90411	L06408-0005
201-90411	L06408-0006
201-90411	L06408-0008
201-90411	L06470-0001
201-90411	L06470-0005
201-90411	L06471-0001
201-90411	L06471-0002
201-90411	L06471-0003
201-90411	L06471-0004
201-90411	L06471-0005
201-90411	L06471-0006
201-90411	L06471-0007
201-90411	L06471-0008
201-90411	L06471-0009
201-90411	L06471-0010
201-90411	L06472-0001
201-90411	L06472-0002
201-90411	L06472-0004
201-90411	L06472-0005
201-90411	L06472-0006
201-90411	L06472-0009
201-90411	L06472-0010
201-90411	L06474-0001
201-90411	L06474-0002
201-90411	L06474-0004
201-90411	L06474-0005
201-90411	L06474-0006
201-90411	L06474-0009
201-90411	L06480-0001

Model Number	Serial Number
201-90411	L06480-0002
201-90411	L06480-0003
201-90411	L06480-0004
201-90411	L06480-0008
201-90411	L06480-0009
201-90411	L06480-0010
201-90411	L06481-0001
201-90411	L06481-0002
201-90411	L06481-0003
201-90411	L06481-0004
201-90411	L06481-0005
201-90411	L06481-0006
201-90411	L06481-0007
201-90411	L06481-0008
201-90411	L06481-0009
201-90411	L06481-0010
201-90411	L06482-0007
201-90411	L06482-0009
201-90411	L06482-0010
201-90411	L06484-0001
201-90411	L06484-0002
201-90411	L06484-0003
201-90411	L06484-0004
201-90411	L06484-0005
201-90411	L06484-0006
201-90411	L06484-0009
201-90411	L06487-0003
201-90411	L06487-0004
201-90411	L06487-0005
201-90411	L06487-0006
201-90411	L06487-0007
201-90411	L06487-0008
201-90411	L06488-0003
201-90411	L06488-0007
201-90411	L06488-0008
201-90411	L06488-0009
201-90411	L06488-0010
201-90411	L06490-0001
201-90411	L06490-0002
201-90411	L06490-0003
201-90411	L06490-0004
201-90411	L06490-0005
201-90411	L06490-0006
201-90411	L06490-0007
201-90411	L06490-0009
201-90411	L06490-0010

Model Number	Serial Number
201-90411	L06491-0001
201-90411	L06491-0002
201-90411	L06491-0003
201-90411	L06491-0005
201-90411	L06491-0006
201-90411	L06491-0007
201-90411	L06493-0002
201-90411	L06493-0004
201-90411	L06493-0010
201-90411	L06494-0007
201-90411	L06496-0001
201-90411	L06496-0002
201-90411	L06496-0003
201-90411	L06496-0004
201-90411	L06496-0005
201-90411	L06496-0006
201-90411	L06496-0007
201-90411	L06496-0008
201-90411	L06496-0009
201-90411	L06496-0010
201-90411	L06497-0001
201-90411	L06497-0002
201-90411	L06497-0003
201-90411	L06497-0004
201-90411	L06497-0008
201-90411	L06497-0009
201-90411	L06521-0002
201-90411	L06521-0003
201-90411	L06521-0004
201-90411	L06521-0006
201-90411	L06521-0007
201-90411	L06521-0008
201-90411	L06521-0009
201-90411	L06522-0009
201-90411	L06523-0008
201-90411	L06523-0010
201-90411	L06524-0001
201-90411	L06524-0002
201-90411	L06524-0003
201-90411	L06524-0004
201-90411	L06524-0005
201-90411	L06524-0006
201-90411	L06524-0007
201-90411	L06524-0008
201-90411	L06524-0009
201-90411	L06524-0010

Appendix A: Specific Model and Serial Numbers in Scope of Notification

CentriMag Console Serial Number Range: L06525-0002 to L06623-0008

Model Number	Serial Number
201-90411	L06525-0002
201-90411	L06525-0003
201-90411	L06525-0006
201-90411	L06527-0001
201-90411	L06527-0002
201-90411	L06527-0003
201-90411	L06527-0005
201-90411	L06527-0007
201-90411	L06527-0010
201-90411	L06528-0001
201-90411	L06529-0001
201-90411	L06529-0003
201-90411	L06529-0004
201-90411	L06529-0005
201-90411	L06529-0006
201-90411	L06529-0007
201-90411	L06529-0008
201-90411	L06529-0009
201-90411	L06529-0010
201-90411	L06553-0002
201-90411	L06556-0008
201-90411	L06557-0001
201-90411	L06558-0008
201-90411	L06581-0001
201-90411	L06581-0002
201-90411	L06581-0003
201-90411	L06581-0004
201-90411	L06581-0005
201-90411	L06581-0006
201-90411	L06581-0007
201-90411	L06581-0008
201-90411	L06582-0003
201-90411	L06582-0004
201-90411	L06582-0005
201-90411	L06582-0006
201-90411	L06582-0007
201-90411	L06583-0010
201-90411	L06584-0003
201-90411	L06585-0001
201-90411	L06585-0002
201-90411	L06585-0003
201-90411	L06585-0004
201-90411	L06585-0006
201-90411	L06585-0007
201-90411	L06585-0008
201-90411	L06587-0002

Model Number	Serial Number
201-90411	L06588-0001
201-90411	L06589-0001
201-90411	L06589-0002
201-90411	L06589-0003
201-90411	L06589-0004
201-90411	L06589-0005
201-90411	L06589-0006
201-90411	L06589-0007
201-90411	L06589-0008
201-90411	L06589-0009
201-90411	L06589-0010
201-90411	L06590-0007
201-90411	L06590-0008
201-90411	L06590-0009
201-90411	L06590-0010
201-90411	L06591-0001
201-90411	L06591-0002
201-90411	L06591-0003
201-90411	L06591-0004
201-90411	L06591-0005
201-90411	L06591-0006
201-90411	L06591-0007
201-90411	L06592-0001
201-90411	L06592-0002
201-90411	L06592-0003
201-90411	L06592-0005
201-90411	L06592-0006
201-90411	L06592-0007
201-90411	L06592-0008
201-90411	L06592-0009
201-90411	L06593-0003
201-90411	L06593-0009
201-90411	L06593-0010
201-90411	L06594-0001
201-90411	L06594-0003
201-90411	L06594-0004
201-90411	L06594-0005
201-90411	L06594-0006
201-90411	L06594-0007
201-90411	L06594-0008
201-90411	L06594-0009
201-90411	L06594-0010
201-90411	L06597-0001
201-90411	L06597-0002
201-90411	L06597-0003
201-90411	L06597-0004

Model Number	Serial Number
201-90411	L06597-0005
201-90411	L06597-0006
201-90411	L06597-0007
201-90411	L06597-0008
201-90411	L06597-0010
201-90411	L06616-0001
201-90411	L06616-0006
201-90411	L06616-0007
201-90411	L06616-0008
201-90411	L06616-0009
201-90411	L06617-0001
201-90411	L06617-0002
201-90411	L06617-0005
201-90411	L06617-0006
201-90411	L06617-0007
201-90411	L06617-0008
201-90411	L06617-0010
201-90411	L06619-0001
201-90411	L06619-0002
201-90411	L06619-0003
201-90411	L06619-0004
201-90411	L06619-0005
201-90411	L06619-0009
201-90411	L06620-0004
201-90411	L06620-0005
201-90411	L06620-0008
201-90411	L06620-0009
201-90411	L06620-0010
201-90411	L06621-0001
201-90411	L06621-0002
201-90411	L06621-0003
201-90411	L06621-0004
201-90411	L06621-0005
201-90411	L06621-0006
201-90411	L06621-0007
201-90411	L06621-0008
201-90411	L06621-0009
201-90411	L06621-0010
201-90411	L06622-0003
201-90411	L06622-0007
201-90411	L06622-0008
201-90411	L06623-0001
201-90411	L06623-0002
201-90411	L06623-0003
201-90411	L06623-0005
201-90411	L06623-0008

Appendix A: Specific Model and Serial Numbers in Scope of Notification

CentriMag Console Serial Number Range: L06623-0009 to L06934-0004

Model Number	Serial Number
201-90411	L06623-0009
201-90411	L06623-0010
201-90411	L06624-0001
201-90411	L06624-0002
201-90411	L06624-0003
201-90411	L06624-0004
201-90411	L06624-0007
201-90411	L06625-0001
201-90411	L06625-0002
201-90411	L06625-0003
201-90411	L06625-0004
201-90411	L06625-0005
201-90411	L06625-0006
201-90411	L06625-0007
201-90411	L06625-0008
201-90411	L06625-0010
201-90411	L06756-0001
201-90411	L06756-0002
201-90411	L06756-0003
201-90411	L06756-0005
201-90411	L06756-0006
201-90411	L06756-0007
201-90411	L06756-0008
201-90411	L06756-0009
201-90411	L06756-0010
201-90411	L06774-0002
201-90411	L06774-0003
201-90411	L06774-0004
201-90411	L06774-0005
201-90411	L06774-0006
201-90411	L06774-0008
201-90411	L06774-0009
201-90411	L06774-0010
201-90411	L06775-0009
201-90411	L06775-0010
201-90411	L06777-0001
201-90411	L06777-0002
201-90411	L06777-0003
201-90411	L06777-0005
201-90411	L06777-0006
201-90411	L06777-0007
201-90411	L06779-0001
201-90411	L06779-0002
201-90411	L06779-0003
201-90411	L06779-0005
201-90411	L06779-0006

Model Number	Serial Number
201-90411	L06779-0008
201-90411	L06779-0009
201-90411	L06779-0010
201-90411	L06781-0001
201-90411	L06781-0002
201-90411	L06781-0003
201-90411	L06800-0001
201-90411	L06800-0003
201-90411	L06800-0005
201-90411	L06800-0007
201-90411	L06800-0008
201-90411	L06800-0009
201-90411	L06800-0010
201-90411	L06801-0008
201-90411	L06801-0009
201-90411	L06802-0007
201-90411	L06803-0001
201-90411	L06803-0003
201-90411	L06803-0004
201-90411	L06803-0006
201-90411	L06803-0007
201-90411	L06803-0008
201-90411	L06803-0009
201-90411	L06805-0009
201-90411	L06806-0001
201-90411	L06806-0002
201-90411	L06806-0003
201-90411	L06806-0004
201-90411	L06806-0005
201-90411	L06806-0006
201-90411	L06806-0007
201-90411	L06806-0008
201-90411	L06806-0009
201-90411	L06806-0010
201-90411	L06807-0001
201-90411	L06807-0002
201-90411	L06807-0003
201-90411	L06807-0004
201-90411	L06807-0005
201-90411	L06807-0008
201-90411	L06807-0009
201-90411	L06807-0010
201-90411	L06863-0005
201-90411	L06864-0001
201-90411	L06864-0002
201-90411	L06864-0003

Model Number	Serial Number
201-90411	L06864-0004
201-90411	L06864-0005
201-90411	L06864-0006
201-90411	L06864-0007
201-90411	L06864-0008
201-90411	L06864-0009
201-90411	L06864-0010
201-90411	L06865-0001
201-90411	L06865-0002
201-90411	L06865-0003
201-90411	L06865-0004
201-90411	L06865-0005
201-90411	L06865-0006
201-90411	L06865-0007
201-90411	L06865-0008
201-90411	L06865-0009
201-90411	L06877-0010
201-90411	L06878-0001
201-90411	L06878-0002
201-90411	L06878-0005
201-90411	L06878-0007
201-90411	L06878-0009
201-90411	L06878-0010
201-90411	L06903-0002
201-90411	L06904-0001
201-90411	L06904-0002
201-90411	L06904-0003
201-90411	L06904-0004
201-90411	L06904-0005
201-90411	L06904-0006
201-90411	L06904-0007
201-90411	L06904-0009
201-90411	L06904-0010
201-90411	L06907-0004
201-90411	L06922-0002
201-90411	L06922-0003
201-90411	L06922-0004
201-90411	L06922-0006
201-90411	L06922-0010
201-90411	L06925-0001
201-90411	L06925-0002
201-90411	L06925-0003
201-90411	L06925-0005
201-90411	L06925-0006
201-90411	L06926-0007
201-90411	L06934-0004

Appendix A: Specific Model and Serial Numbers in Scope of Notification

CentriMag Console Serial Number Range: L06934-0005 to L07459-0001

Model Number	Serial Number
201-90411	L06934-0005
201-90411	L06934-0006
201-90411	L06934-0008
201-90411	L06934-0010
201-90411	L06940-0003
201-90411	L06941-0005
201-90411	L06946-0010
201-90411	L06947-0002
201-90411	L06947-0003
201-90411	L06947-0004
201-90411	L06947-0006
201-90411	L06947-0007
201-90411	L06947-0008
201-90411	L06949-0002
201-90411	L06949-0005
201-90411	L06949-0007
201-90411	L07020-0001
201-90411	L07021-0001
201-90411	L07026-0001
201-90411	L07032-0001
201-90411	L07033-0001
201-90411	L07034-0001
201-90411	L07038-0001
201-90411	L07040-0001
201-90411	L07041-0001
201-90411	L07042-0001
201-90411	L07045-0001
201-90411	L07046-0001
201-90411	L07047-0001
201-90411	L07077-0001
201-90411	L07078-0001
201-90411	L07080-0001
201-90411	L07081-0001
201-90411	L07131-0001
201-90411	L07132-0001
201-90411	L07133-0001
201-90411	L07134-0001
201-90411	L07135-0001
201-90411	L07136-0001
201-90411	L07137-0001
201-90411	L07138-0001
201-90411	L07139-0001
201-90411	L07140-0001
201-90411	L07143-0001
201-90411	L07145-0001
201-90411	L07146-0001

Model Number	Serial Number
201-90411	L07147-0001
201-90411	L07148-0001
201-90411	L07152-0001
201-90411	L07154-0001
201-90411	L07159-0001
201-90411	L07160-0001
201-90411	L07165-0001
201-90411	L07166-0001
201-90411	L07167-0001
201-90411	L07168-0001
201-90411	L07169-0001
201-90411	L07170-0001
201-90411	L07171-0001
201-90411	L07172-0001
201-90411	L07175-0001
201-90411	L07176-0001
201-90411	L07178-0001
201-90411	L07179-0001
201-90411	L07180-0001
201-90411	L07181-0001
201-90411	L07183-0001
201-90411	L07184-0001
201-90411	L07186-0001
201-90411	L07187-0001
201-90411	L07190-0001
201-90411	L07191-0001
201-90411	L07192-0001
201-90411	L07193-0001
201-90411	L07194-0001
201-90411	L07195-0001
201-90411	L07196-0001
201-90411	L07197-0001
201-90411	L07199-0001
201-90411	L07200-0001
201-90411	L07201-0001
201-90411	L07203-0001
201-90411	L07204-0001
201-90411	L07206-0001
201-90411	L07207-0001
201-90411	L07208-0001
201-90411	L07209-0001
201-90411	L07234-0001
201-90411	L07238-0001
201-90411	L07247-0001
201-90411	L07251-0001
201-90411	L07257-0001

Model Number	Serial Number
201-90411	L07264-0001
201-90411	L07268-0001
201-90411	L07269-0001
201-90411	L07282-0001
201-90411	L07285-0001
201-90411	L07291-0001
201-90411	L07298-0001
201-90411	L07299-0001
201-90411	L07300-0001
201-90411	L07303-0001
201-90411	L07306-0001
201-90411	L07310-0001
201-90411	L07324-0001
201-90411	L07328-0001
201-90411	L07346-0001
201-90411	L07360-0001
201-90411	L07377-0001
201-90411	L07385-0001
201-90411	L07390-0001
201-90411	L07393-0001
201-90411	L07402-0001
201-90411	L07405-0001
201-90411	L07406-0001
201-90411	L07412-0001
201-90411	L07433-0001
201-90411	L07437-0001
201-90411	L07438-0001
201-90411	L07439-0001
201-90411	L07440-0001
201-90411	L07442-0001
201-90411	L07443-0001
201-90411	L07444-0001
201-90411	L07445-0001
201-90411	L07446-0001
201-90411	L07447-0001
201-90411	L07448-0001
201-90411	L07449-0001
201-90411	L07450-0001
201-90411	L07451-0001
201-90411	L07452-0001
201-90411	L07453-0001
201-90411	L07454-0001
201-90411	L07456-0001
201-90411	L07457-0001
201-90411	L07458-0001
201-90411	L07459-0001

Appendix A: Specific Model and Serial Numbers in Scope of Notification

CentriMag Console Serial Number Range: L07460-0001 to L07874-0001

Model Number	Serial Number
201-90411	L07460-0001
201-90411	L07461-0001
201-90411	L07465-0001
201-90411	L07466-0001
201-90411	L07467-0001
201-90411	L07468-0001
201-90411	L07472-0001
201-90411	L07473-0001
201-90411	L07475-0001
201-90411	L07477-0001
201-90411	L07478-0001
201-90411	L07479-0001
201-90411	L07495-0001
201-90411	L07538-0001
201-90411	L07552-0001
201-90411	L07553-0001
201-90411	L07562-0001
201-90411	L07565-0001
201-90411	L07568-0001
201-90411	L07571-0001
201-90411	L07574-0001
201-90411	L07575-0001
201-90411	L07576-0001
201-90411	L07590-0001
201-90411	L07593-0001
201-90411	L07602-0001
201-90411	L07603-0001
201-90411	L07605-0001
201-90411	L07608-0001
201-90411	L07613-0001
201-90411	L07615-0001
201-90411	L07616-0001
201-90411	L07617-0001
201-90411	L07618-0001
201-90411	L07619-0001
201-90411	L07622-0001
201-90411	L07623-0001
201-90411	L07657-0001
201-90411	L07659-0001
201-90411	L07660-0001
201-90411	L07662-0001
201-90411	L07663-0001
201-90411	L07664-0001
201-90411	L07665-0001
201-90411	L07666-0001
201-90411	L07667-0001

Model Number	Serial Number
201-90411	L07669-0001
201-90411	L07672-0001
201-90411	L07674-0001
201-90411	L07675-0001
201-90411	L07677-0001
201-90411	L07678-0001
201-90411	L07679-0001
201-90411	L07680-0001
201-90411	L07681-0001
201-90411	L07682-0001
201-90411	L07683-0001
201-90411	L07684-0001
201-90411	L07685-0001
201-90411	L07687-0001
201-90411	L07688-0001
201-90411	L07689-0001
201-90411	L07690-0001
201-90411	L07691-0001
201-90411	L07692-0001
201-90411	L07693-0001
201-90411	L07694-0001
201-90411	L07695-0001
201-90411	L07696-0001
201-90411	L07697-0001
201-90411	L07698-0001
201-90411	L07699-0001
201-90411	L07700-0001
201-90411	L07771-0001
201-90411	L07772-0001
201-90411	L07773-0001
201-90411	L07774-0001
201-90411	L07775-0001
201-90411	L07776-0001
201-90411	L07777-0001
201-90411	L07778-0001
201-90411	L07779-0001
201-90411	L07780-0001
201-90411	L07781-0001
201-90411	L07782-0001
201-90411	L07783-0001
201-90411	L07785-0001
201-90411	L07787-0001
201-90411	L07788-0001
201-90411	L07789-0001
201-90411	L07790-0001
201-90411	L07791-0001

Model Number	Serial Number
201-90411	L07792-0001
201-90411	L07793-0001
201-90411	L07794-0001
201-90411	L07795-0001
201-90411	L07796-0001
201-90411	L07797-0001
201-90411	L07798-0001
201-90411	L07799-0001
201-90411	L07800-0001
201-90411	L07803-0001
201-90411	L07804-0001
201-90411	L07805-0001
201-90411	L07806-0001
201-90411	L07807-0001
201-90411	L07810-0001
201-90411	L07811-0001
201-90411	L07816-0001
201-90411	L07817-0001
201-90411	L07819-0001
201-90411	L07820-0001
201-90411	L07821-0001
201-90411	L07823-0001
201-90411	L07824-0001
201-90411	L07825-0001
201-90411	L07826-0001
201-90411	L07827-0001
201-90411	L07828-0001
201-90411	L07829-0001
201-90411	L07830-0001
201-90411	L07831-0001
201-90411	L07832-0001
201-90411	L07833-0001
201-90411	L07834-0001
201-90411	L07835-0001
201-90411	L07846-0001
201-90411	L07847-0001
201-90411	L07852-0001
201-90411	L07853-0001
201-90411	L07861-0001
201-90411	L07864-0001
201-90411	L07867-0001
201-90411	L07869-0001
201-90411	L07871-0001
201-90411	L07872-0001
201-90411	L07873-0001
201-90411	L07874-0001

Appendix A: Specific Model and Serial Numbers in Scope of Notification

CentriMag Console Serial Number Range: L07876-0001 to L08002-0001,
19258710 to 19605455, L01627-0007 to L05212-0001

Model Number	Serial Number
201-90411	L07876-0001
201-90411	L07877-0001
201-90411	L07878-0001
201-90411	L07881-0001
201-90411	L07887-0001
201-90411	L07889-0001
201-90411	L07891-0001
201-90411	L07893-0001
201-90411	L07894-0001
201-90411	L07895-0001
201-90411	L07896-0001
201-90411	L07897-0001
201-90411	L07898-0001
201-90411	L07899-0001
201-90411	L07900-0001
201-90411	L07902-0001
201-90411	L07903-0001
201-90411	L07904-0001
201-90411	L07905-0001
201-90411	L07908-0001
201-90411	L07910-0001
201-90411	L07911-0001
201-90411	L07912-0001
201-90411	L07913-0001
201-90411	L07914-0001
201-90411	L07915-0001
201-90411	L07916-0001
201-90411	L07919-0001
201-90411	L07920-0001
201-90411	L07932-0001
201-90411	L07934-0001
201-90411	L07935-0001
201-90411	L07938-0001
201-90411	L07940-0001
201-90411	L07941-0001
201-90411	L07942-0001
201-90411	L07943-0001
201-90411	L07945-0001
201-90411	L07947-0001
201-90411	L07950-0001
201-90411	L07951-0001
201-90411	L07952-0001
201-90411	L07953-0001
201-90411	L07954-0001
201-90411	L07955-0001
201-90411	L07958-0001

Model Number	Serial Number
201-90411	L07961-0001
201-90411	L07962-0001
201-90411	L07963-0001
201-90411	L07964-0001
201-90411	L07966-0001
201-90411	L07968-0001
201-90411	L07969-0001
201-90411	L07971-0001
201-90411	L07973-0001
201-90411	L07974-0001
201-90411	L07976-0001
201-90411	L07977-0001
201-90411	L07978-0001
201-90411	L07979-0001
201-90411	L07981-0001
201-90411	L07982-0001
201-90411	L07983-0001
201-90411	L07984-0001
201-90411	L07987-0001
201-90411	L07989-0001
201-90411	L07990-0001
201-90411	L07991-0001
201-90411	L07994-0001
201-90411	L07995-0001
201-90411	L07996-0001
201-90411	L07997-0001
201-90411	L07999-0001
201-90411	L08000-0001
201-90411	L08002-0001
L201-90411	19258710
L201-90411	19280695
L201-90411	19280702
L201-90411	19476967
L201-90411	19604016
L201-90411	19605455
L201-90411	L01627-0007
L201-90411	L01627-0009
L201-90411	L01965-0001
L201-90411	L01965-0009
L201-90411	L02293-0002
L201-90411	L02305-0003
L201-90411	L03009-0006
L201-90411	L03009-0007
L201-90411	L03023-0008
L201-90411	L03065-0008
L201-90411	L03075-0005

Model Number	Serial Number
L201-90411	L03075-0007
L201-90411	L03075-0008
L201-90411	L03075-0010
L201-90411	L03082-0001
L201-90411	L03117-0004
L201-90411	L03117-0007
L201-90411	L03117-0010
L201-90411	L03123-0008
L201-90411	L03226-0001
L201-90411	L03226-0002
L201-90411	L03226-0004
L201-90411	L03314-0003
L201-90411	L03530-0004
L201-90411	L03530-0008
L201-90411	L03543-0002
L201-90411	L03659-0008
L201-90411	L03659-0010
L201-90411	L03812-0008
L201-90411	L03812-0009
L201-90411	L03849-0001
L201-90411	L03849-0002
L201-90411	L03849-0005
L201-90411	L03882-0006
L201-90411	L03882-0009
L201-90411	L03897-0003
L201-90411	L03897-0004
L201-90411	L03897-0006
L201-90411	L04503-0003
L201-90411	L04503-0004
L201-90411	L04503-0005
L201-90411	L04696-0005
L201-90411	L05083-0003
L201-90411	L05083-0004
L201-90411	L05083-0005
L201-90411	L05083-0006
L201-90411	L05083-0007
L201-90411	L05096-0005
L201-90411	L05116-0008
L201-90411	L05150-0004
L201-90411	L05150-0006
L201-90411	L05150-0008
L201-90411	L05150-0009
L201-90411	L05168-0007
L201-90411	L05168-0008
L201-90411	L05202-0008
L201-90411	L05212-0001

Appendix A: Specific Model and Serial Numbers in Scope of Notification

CentriMag Console Serial Number Range: L05229-0007 to L07980-0001

Model Number	Serial Number
L201-90411	L05229-0007
L201-90411	L05238-0001
L201-90411	L05238-0002
L201-90411	L05238-0003
L201-90411	L05569-0007
L201-90411	L05569-0010
L201-90411	L05582-0006
L201-90411	L05582-0009
L201-90411	L05618-0002
L201-90411	L05724-0003
L201-90411	L05724-0008
L201-90411	L05729-0001
L201-90411	L05837-0009
L201-90411	L05851-0001
L201-90411	L05869-0008
L201-90411	L05869-0009
L201-90411	L05869-0010
L201-90411	L05877-0001
L201-90411	L05877-0009
L201-90411	L05877-0010
L201-90411	L05917-0008
L201-90411	L05963-0006
L201-90411	L05963-0007
L201-90411	L05963-0008
L201-90411	L06302-0008
L201-90411	L06338-0007
L201-90411	L06338-0009
L201-90411	L06355-0001
L201-90411	L06355-0003
L201-90411	L06359-0004
L201-90411	L06364-0004
L201-90411	L06364-0006
L201-90411	L06379-0001
L201-90411	L06379-0002
L201-90411	L06379-0003
L201-90411	L06379-0004
L201-90411	L06379-0006
L201-90411	L06380-0002
L201-90411	L06381-0005
L201-90411	L06381-0006
L201-90411	L06381-0008
L201-90411	L06381-0009
L201-90411	L06385-0006
L201-90411	L06385-0007
L201-90411	L06385-0008
L201-90411	L06386-0003

Model Number	Serial Number
L201-90411	L06386-0010
L201-90411	L06387-0005
L201-90411	L06387-0006
L201-90411	L06388-0001
L201-90411	L06388-0003
L201-90411	L06388-0004
L201-90411	L06401-0002
L201-90411	L06401-0003
L201-90411	L06401-0005
L201-90411	L06403-0001
L201-90411	L06406-0009
L201-90411	L06470-0003
L201-90411	L06470-0004
L201-90411	L06474-0007
L201-90411	L06474-0010
L201-90411	L06483-0002
L201-90411	L06483-0004
L201-90411	L06483-0007
L201-90411	L06484-0007
L201-90411	L06484-0008
L201-90411	L06487-0001
L201-90411	L06488-0001
L201-90411	L06488-0002
L201-90411	L06488-0004
L201-90411	L06493-0001
L201-90411	L06493-0003
L201-90411	L06493-0006
L201-90411	L06493-0008
L201-90411	L06497-0005
L201-90411	L06497-0006
L201-90411	L06497-0010
L201-90411	L06521-0005
L201-90411	L06522-0010
L201-90411	L06523-0004
L201-90411	L06523-0009
L201-90411	L06527-0008
L201-90411	L06527-0009
L201-90411	L06557-0005
L201-90411	L06585-0010
L201-90411	L06616-0003
L201-90411	L06616-0004
L201-90411	L06616-0010
L201-90411	L06619-0008
L201-90411	L06622-0006
L201-90411	L06622-0009
L201-90411	L06624-0005

Model Number	Serial Number
L201-90411	L06775-0004
L201-90411	L06775-0005
L201-90411	L06777-0008
L201-90411	L06779-0004
L201-90411	L06800-0002
L201-90411	L06802-0010
L201-90411	L06807-0006
L201-90411	L06865-0010
L201-90411	L06877-0007
L201-90411	L06877-0008
L201-90411	L06877-0009
L201-90411	L06925-0010
L201-90411	L06934-0009
L201-90411	L06941-0001
L201-90411	L06941-0002
L201-90411	L06941-0003
L201-90411	L06941-0004
L201-90411	L07025-0001
L201-90411	L07286-0001
L201-90411	L07287-0001
L201-90411	L07293-0001
L201-90411	L07349-0001
L201-90411	L07357-0001
L201-90411	L07476-0001
L201-90411	L07489-0001
L201-90411	L07661-0001
L201-90411	L07784-0001
L201-90411	L07814-0001
L201-90411	L07815-0001
L201-90411	L07842-0001
L201-90411	L07866-0001
L201-90411	L07875-0001
L201-90411	L07901-0001
L201-90411	L07907-0001
L201-90411	L07980-0001